


# iRESPITRONICS

Advancing respiratory care *through technology*



## Elite Plus

Where Boundless Innovation Meets the Symphony of Clear Breathing:  
Elite Econo Unleashes Unrivaled High-Tech Respiratory Solutions.

# Navigating the Features



## Single Controller

The iRespitronics Elite Plus is powered by a single controller which uses either an A/C adaptor or 2x AA batteries.



## Heat-Free Operation

iRespitronics employs advanced technology that ensures heat is not produced during operation which maintains the integrity and effectiveness of the medication.



## Cost Saver

The Elite Plus with the Controller and flow sensor senses the patient's breathing pattern and automatically initiates nebulization ONLY when the patient inhales.



## Connection to Face Mask

The iRespitronics Elite Plus can connect directly to facemasks using standard connections common at most hospitals.



## Connect to Tracheotomy Tube

The iRespitronics Elite Plus can connect directly to a tracheotomy tube using standard connections common at most hospitals.



## Connect to Endotracheal Tube

The iRespitronics Elite Plus can connect directly to an endotracheal tube using standard connections common at most hospitals.



## Nebulization Rate

Nebulization rate refers to the speed at which the reservoir transforms liquid into mist.  
>12 mL/Hour



## Vibrating Mesh

200 Particle Size Per Minute  
2mL Nebulizer Rate



## Auto Shut-Off

The Elite Plus and Controller will turn off automatically when the medication in the reservoir is out.



## LED Color-State

The Elite Plus when connected to the Controller signals the level of medication in the reservoir with the colors Green, Yellow, and Red.

# Navigating the Functionality

## Reservoir

The medication is filled into the reservoir. The reservoir's capacity is 15 mL.



## Ventilator Tubing

The Elite Plus seamlessly connects to standard hospital ventilator tubing.

## Air Input

The Air Input is where the Ventilator Tubing pushing air as the patient inhales

## Electrical Connection

Connecting electrical wire from Controller to the Respitronics Elite Plus to powerup the device.

# Elite Pro: Enhanced with Integrated Controller and Flow Sensor for the Ultimate Respiratory Powerhouse



## A/C, D/C, & USB Powered

The iRespironics Pro is powered by a single controller which offers versatile power options: A/C adapter, USB adapter or 2x AA Batteries.



## LED Color-State

The Elite Pro when connected to the Controller signals the level of medication in the reservoir with the colors Green, Yellow, and Red.



## Cost Saver

The Elite Plus with the Controller and flow sensor senses the patient's breathing pattern and automatically initiates nebulization ONLY when the patient inhales.



## Auto Shut-Off

The Elite Pro and Controller will turn off automatically when the medication in the reservoir is out.